



THE UNIVERSITY OF ALABAMA

NAA techniques and applications

Ryan MacLellan

LRT 2010

August 28, 2010



Techniques typically used to characterize natural radioactivity

Direct γ -ray counting

- Sensitive to the background creating isotopes of the decay series directly.
- Requires long counting time and large sample mass to achieve good sensitivity.

10^{-14} g/g U produces only 1 decay per day per 100kg of sample mass.

Mass spectroscopy

- Sensitivity limits of:
 - 10^{-9} g/g K
 - $10^{-12/-13}$ g/g U and Th
- Requires sample to be digested by acids (ICPMS) or conductors (GDMS).




Neutron Activation Analysis

- Sensitivity limits of:
 - 10^{-9} g/g K
 - $10^{-12/-13}$ g/g U and Th
- High sensitivity measurement can take up to one month due to long half life of Th activation product: ^{233}Pa .
- Only applicable to matrix materials with low neutron capture cross sections (most metals not suitable).




Gold plated silicon wafer (etched for ICP-MS)

2d after neutron activation



1 month after neutron activation



$$t_{1/2} = 2.7\text{d}$$



Neutron Activation Analysis

- Sensitivity can be increased by removing uninteresting side activities before counting.
- For liquid scintillator samples, post-radiation chemistry increased sensitivity to 10^{-15} g/g U/Th at UA.
- A group at TU Munich achieved 10^{-16} g/g U/Th utilizing chemistry and an onsite reactor.
- Only feasible for screening limited number of materials as each requires custom tailored chemical processes.



Neutron Activation Analysis

- Any neutron activation analysis requires a sample handling facility to handle considerable amounts of open radioactivity.
- Consistent pre-analysis treatment is critical to avoid surface contamination
 - At UA this is done in a class 500 clean room
 - High purity solvents, acids, and water for cleaning
- Activation costs for UA 1000--2000USD.
- Sample counting on, ideally multiple, large Ge detectors to combat decaying signal. Low background detectors are not required (except RNAA).



Procedure at UA

- Sample preparation in class 500 clean room.
- All lab supplies etched with ultra pure acid and rinsed with purified water.
- Depending on the required irradiation duration, sample activated in a small PE bottle or miniature quartz vial.





A few examples of reactors


Thermal neutron flux of some reactors with pneumatic sample insertion:

Reactor	Location	Thermal n Flux ($\text{cm}^{-2}\text{s}^{-1}$)
MITR	MIT	5×10^{13}
HFIR	Oak Ridge	4×10^{14}
MURR	University of Missouri	6×10^{14}
RFM II	TU Munich	8×10^{14}



August 28, 2010

Ryan MacLellan, LRT 2010






August 28, 2010

Ryan MacLellan, LRT 2010

10






NA Analysis








Evidence of natural radioactivity tracers





Time evolution of measured activity





Simply calculating constituents

Number of activated nuclei

$$A = \underbrace{\phi_n \sigma_{\text{tot}} N (1 - e^{-\lambda t_{\text{irr}}})}_{\text{Number of activated nuclei}} \times \underbrace{e^{-\lambda t}}_{\text{Fraction remaining at time } t}$$

Concentration of parent nuclei: $X[\text{ng/g}] = 10^9 \frac{N}{N_A} \frac{M}{x} / m_{\text{sample}}$



Simply calculating constituents

Number of activated nuclei

$$A = \underbrace{\phi_n \sigma_{\text{tot}} N (1 - e^{-\lambda t_{\text{irr}}})}_{\text{Number of activated nuclei}} \times \underbrace{e^{-\lambda t}}_{\text{Fraction remaining at time } t}$$

Concentration of parent nuclei: $X[\text{ng/g}] = 10^9 \frac{N}{N_A} \frac{M}{x} / m_{\text{sample}}$

For a two stage decay things get more complicated but still straight forwardly:

$$A = \phi_n \sigma_{\text{tot}} N \left(\frac{\lambda'}{\lambda' - \lambda} \frac{1 - e^{-\lambda t_{\text{irr}}}}{\lambda} e^{-\lambda t} - \frac{\lambda}{\lambda' - \lambda} \frac{1 - e^{-\lambda' t_{\text{irr}}}}{\lambda'} e^{-\lambda' t} \right) \lambda$$



Simply calculating constituents

Number of activated nuclei

$$A = \underbrace{\phi_n \sigma_{\text{tot}} N (1 - e^{-\lambda t_{\text{irr}}})}_{\text{Number of activated nuclei}} \times \underbrace{e^{-\lambda t}}_{\text{Fraction remaining at time } t}$$

Concentration of parent nuclei: $X[\text{ng/g}] = 10^9 \frac{N}{N_A} \frac{M}{x} / m_{\text{sample}}$

For a two stage decay things get more complicated but still straight forwardly:

$$A = \phi_n \sigma_{\text{tot}} N \left(\frac{\lambda'}{\lambda' - \lambda} \frac{1 - e^{-\lambda t_{\text{irr}}}}{\lambda} e^{-\lambda t} - \frac{\lambda}{\lambda' - \lambda} \frac{1 - e^{-\lambda' t_{\text{irr}}}}{\lambda'} e^{-\lambda' t} \right) \lambda$$



Simply calculating constituents

Number of activated nuclei

$$A = \underbrace{\phi_n \sigma_{\text{tot}} N (1 - e^{-\lambda t_{\text{irr}}})}_{\text{Number of activated nuclei}} \times \underbrace{e^{-\lambda t}}_{\text{Fraction remaining at time } t}$$


Concentration of parent nuclei: $X[\text{ng/g}] = 10^9 \frac{N}{N_A} \frac{M}{x} / m_{\text{sample}}$

For a two stage decay things get more complicated but still straight forwardly:


$$A = \phi_n \sigma_{\text{tot}} N \left(\frac{\lambda'}{\lambda' - \lambda} \frac{1 - e^{-\lambda t_{\text{irr}}}}{\lambda} e^{-\lambda t} - \frac{\lambda}{\lambda' - \lambda} \frac{1 - e^{-\lambda' t_{\text{irr}}}}{\lambda'} e^{-\lambda' t} \right) \lambda$$




Parametrize the neutron flux: thermal vs. epithermal neutrons



At
MITR and HFIR
 $T_m \approx 340\text{K}$



Parametrize the product of neutron flux and cross section; thermal vs. epithermal neutrons





Parametrize the product of neutron flux and cross section; thermal vs. epithermal neutrons


$$\int \sigma(E)\Phi(E)dE \approx \phi_{\text{th}} \left(\langle \sigma_{\text{th}} \rangle + \frac{1}{f} \langle \sigma_{\text{epi}} \rangle \right) \equiv \phi_{\text{n}} \sigma_{\text{tot}}$$

Thermal neutron flux: $\phi_{\text{n}} \equiv \phi_{\text{th}}$

Epithermal flux suppression factor: $f \equiv \frac{\phi_{\text{th}}}{\phi_{\text{epi}}}$

$$\langle \sigma_{\text{th/epi}} \rangle = \frac{\int_{E_{\text{th/epi}}} \sigma(E)\Phi(E)dE}{\int_{E_{\text{th/epi}}} \Phi(E)dE}$$

$$\phi_{\text{th/epi}} = \int_{E_{\text{th/epi}}} \Phi(E)dE$$





National Institute of Standards & Technology

Certificate of Analysis

Standard Reference Material[®] 1633b

Constituent Elements in Coal Fly Ash

Al	15.05	±	0.27
Ca	1.51	±	0.06
Fe	7.78	±	0.23
Mg	0.482	±	0.008
K	1.95	±	0.03
Si	23.02	±	0.08
Na	0.201	±	0.003
S	0.2075	±	0.0011
Ti	0.791	±	0.014

As	136.2	±	2.6
Ba	709	±	27
Cd	0.784	±	0.006
Cr	198.2	±	4.7
Cu	112.8	±	2.6
Pb	68.2	±	1.1
Mn	131.8	±	1.7
Hg	0.1431	±	0.0018
Ni	120.6	±	1.8
Se	10.26	±	0.17
Sr	1041	±	14
Th	25.7	±	1.3
U	8.79	±	0.36
V	295.7	±	3.6



Global χ^2 minimization of

$$A = \phi_n \sigma_{\text{tot}} N (1 - e^{-\lambda t_{\text{irr}}}) \times e^{-\lambda t}$$


N measured vs N expected

$$N = N_A \frac{x}{M} m_{\text{sample}} m_{\text{frac}}$$

yields ϕ_n and f




MITR ϕ_n as a function of reactor power





MITR f as a function of reactor power





Summary

- NAA has been successfully used in the preparation of previous successful low background experiments. It offers ppt to sub ppt-sensitivity depending effort made to avoid source related background.
- Complements ICP-MS technique. Well suited for most plastics that are inaccessible to traditional MS.

One final note: document surveys


Systematic study of trace radioactive
impurities in candidate construction materials
for EXO-200


10 pages of table (small font) listing all materials characterization for EXO-200 construction (useful and otherwise). Most cited EXO paper to date!

D. S. Leonard *et al.*, Systematic study of trace radioactive impurities in candidate construction materials for EXO-200, Nucl. Instrum. Meth. A **591**, 490 (2008) [arXiv:0709.4524]



UA KamLAND scintillator analysis





Status	Radionuclide	Activity	Activity error	Model prediction	χ^2 contribution
on	^{24}Na	13640	200	13537.6	0.2
on	^{42}K	6320	110	6508.4	2.8
on	^{59}Fe	2220	130	1810.4	9.7
on	^{51}Cr	1650	60	1601.9	0.6
off	^{76}As	10760	210	13052.4	119.0
on	^{203}Hg	10.5	3.6	0.4	7.9
on	^{233}Pa	654	36	639.4	0.2
on	^{239}Np	770	50	769.7	0.0
off	^{124}Sb	36.2	2.4	44.8	12.3
off	^{46}Sc	5750	60	5715.3	0.4
off	^{60}Co	335.8	1.8	316.9	106.6
off	^{65}Zn	129.9	3.3	93.3	123.0
off	^{82}Br	142	6	121.6	13.3
off	^{140}La	11572	31	12324.2	598.6
off	^{181}Hf	93	10	77.5	2.4
off	^{182}Ta	59.6	2.9	58.5	0.2

Effective Neutron Capture Cross Sections

Andreas Piepke
Department of Physics and Astronomy
University of Alabama

April 18, 2001

Isotope	Tabulated [barn]			This evaluation [barn]		
	σ	$\langle\sigma_{th}\rangle$	$\langle\sigma_{epi}\rangle$	$\langle\sigma_{th}\rangle$	$\langle\sigma_{th}\rangle$	$\langle\sigma_{epi}\rangle$
	293.6 K	293.6 K		293.6 K	340 K	
^{23}Na	0.531	0.471	0.311	0.467	0.434	0.349
^{41}K	1.459	1.294	1.58	1.298	1.206	1.664
^{50}Cr	15.92	14.12	7.38	14.13	13.13	8.30
^{58}Fe	1.300	1.153	1.357	1.020	0.948	1.561
^{59}Co	37.18	32.96	75.51	33.07	30.72	76.68
^{64}Zn				0.677	0.629	1.564
^{81}Br	2.690	2.385	46.63	2.385	2.217	51.03
^{197}Au				88.02	81.99	1566.7
^{232}Th	7.400	6.532	84.35	6.567	6.095	86.91
^{238}U	2.717	2.414	278.1	2.421	2.251	278.3