Ultra-Low Background Counting and Assay Studies At SNOLAB

lan Lawson SNOLAB

2015 CAP Congress University of Alberta







AB

SN



Outline

Motivation for Low Background Counters
Advantages of being deep

Current Facilities in Operation at SNOLAB

- → PGT Ge detector
- Canberra Well detector
- Electrostatic Counters (ESCs)
- → Alpha/Beta Counters
- → Radon Emanation (see talk by C. Jillings)
- → X-ray Fluorescence
- Low Background Database

•Future Low Background Counters and Facilities

Motivation

- Many of the experiments currently searching for dark matter, studying properties of neutrinos or searching for neutrinoless double-beta decay require very low levels of radioactive backgrounds both in their own construction materials and in the surrounding environment.
- These low background levels are required so that the experiments can achieve the required sensitivities for their searches.
- SNOLAB has several facilities which are used to directly measure these radioactive backgrounds.
- The backgrounds in question are on the order of 1 mBq or 1 ppb for ²³⁸U, ²³²Th and ²³⁵U and 1 ppm for ⁴⁰K, or better, measurements down to 1 ppt are required for many components.
- The problem backgrounds can include gammas, alphas and neutrons or resulting interaction products.
- The goal is to measure these backgrounds and then to reduce them to be as low as reasonably achievable.

17/06/2015



Uranium Decay Chain

Uranium – Radium Gamma Intensities			A = 4	hn + 2				63.29 4.84 92.38 2.81 92.80 2.77 112.81 0.28	Th 234 24.10 d	49.55 0.064 113.5 0.010	U 238 4.468x10 ⁹ a		
										1001.03 0.837 766.38 0.294	Pa 234 [*] 1.17 m ^{6.7 h}	2.269 98.2%	
	351.932 37.6 295.224 19.3 241.997 7.43 53.2275 1.2 785.96 1.07	Pb 214 26.8(9) m	α none β none	Po 218 3.10(1) m 9.980% 0.020%	◄ 511 0.076	Rn 222 3.8235(3) d	4 186.211 3.59	Ra 226 1600(1) a	67.672 0.378	Th 230 7.538x10 ⁴ a	53.20 0.123	U 234 7.455x10 ⁵ a	
799 99 298 79 1316 21 1210 17 1070 12 1110 6.9 2010 6.9	Tl 210 1.30(3) m	3 609.312 46.1 3 1764.494 15.4 6 1120.287 15.4 3 1238.110 5.79 3 2204.21 5.08 3 768.356 4.94 3 1377.669 4.00 8 934.061 3.03	α none Bi 214 19.9(4) m 0.276% 99.724%	none	At 218 1.5 s								
	46.539 4.25	Pb 210 22.3(2) a	∢ 799.7 0.0104	Po 214 164.3(20) us									
		none	Bi 210 5.013 d										
		Pb 206 stable	4 803.10 0.00121	Po 210 138.376 d									

17/06/2015

CAP 2015

- 4

Thorium Decay Chain

Thorium A = 4n Gamma Intensities				13.52 1.600 16.2 0.72 12.75 0.304 15.5 0.16	Ra 228 5.75 a	6 3.823 0.264 204.68 0.021	Th 232 1.405x10 ¹⁰ a					
								911.204 25.8 968.971 15.8 338.320 11.27 964.766 4.99 463.004 4.40 794.947 4.25 209.253 3.89	Ac 228 6.15 h			
	238.632 43.3 300.087 3.28 115.183 0.592	Pb 212 10.64(1) h	804.9 0.0019	Po 216 145(2) ms	▼ 549.76 0.114	Rn 220 55.6(1) s	▲ 240.986 4.10	Ra 224 3.66(4) d	84.373 1.220 215.983 0.254 131.613 0.131 166.410 0.104	Th 228 1.9116(16) a		
2614.533 99.0 583.191 84.5 510.77 22.6 860.564 12.42 277.351 6.31 763.13 1.81	Tl 208 3.053(4) m	$\frac{\alpha}{39.858}$ 1.091	Bi 212 60.55(6) m 35.94% 64.06%	β 727.330 6.58 1620.50 1.49 785.37 1.102								
		Pb 208 stable	◄	Po 212 299(2) ns								

SNCLAB

Other Interesting Isotopes

SN



SNOLAB PGT HPGe Counter (The workhorse detector at SNOLAB)



17/06/2015

SN



SNOLAB PGT HPGe Detector Specifications

Motivation

• Survey materials for new, existing and proposed experiments (to be) located @ SNOLAB, such as SNO/SNO+, DEAP/CLEAN, PICASSO/COUPP/PICO, EXO, ... Also survey materials for the DM-ICE and DRIFT experiments, and Canberra.

•Constructed @ SNOLAB in 2005, detector was in UG storage from 1997, continuous operations since 2005

Counter manufactured by PGT in 1992

• Endcap diameter: 83 mm,

- Crystal volume: 210 cm³
- Relative Efficiency is 55% wrt a 7.62 cm dia x 7.62 cm NaI(Tl) detector,
- Resolution 1.8 keV FWHM.

Shielding

- 2 inches Cu + 8 inches Pb
- Nitrogen purge at 2L/min to keep radon out, as the lab radon levels are 150 Bq/m³.

Detection Region

• Energy: 90 – 3000 keV

17/06/2015



PGT HPGe Typical Detector Sensitivity (for a standard 1L or 1 kg sample counted for one week)

Isotope	Sensitivity for Standard Size Samples	Sensitivity for Standard Size Samples
²³⁸ U	0.15 mBq/kg	12 ppt
²³⁵ U	0.15 mBq/kg	264 ppt
²³² Th	0.13 mBq/kg	32 ppt
⁴⁰ K	1.70 mBq/kg	54 ppt
⁶⁰ Co	0.06 mBq/kg	
¹³⁷ Cs	0.17 mBq/kg	
⁵⁴ Mn	0.06 mBq/kg	

SN

AB

SNOLAB

17/06/2015

Measurements To Date For Each Experiment

Experiment or Laboratory	Total (2005 - Today)
SNO	11
SNO+	116
SNOLAB	73
EXO	12
MiniCLEAN	56
DEAP	101
HALO	7
PICASSO	9
DM-ICE / DRIFT	23
COUPP / PICO	61
DAMIC	12
Total	481
Calibrations & Tests	118

Counting time per sample averages 6 days.

Non-experimental measurements include water/snow samples from Fukushima fallout. See publication: Can. J. Phys. 90: 599–603 (2012) CAP 2015





SN



17/06/2015





SNOLAB Canberra Well Detector Specifications

Motivation

• Survey very small quantities of materials, concentrated samples or very expensive materials. Used by DAMIC, DEAP, PICO & SNO+ so far.

•Constructed by Canberra using low activity materials and shielding.

- Counter manufactured by Canberra in 2011 and refurbished in 2012, the cold finger was lengthened as it was too short to fit the shielding and the tail end and crystal holder were replaced to reduce radioactivity levels.
- Crystal volume: 300 cm³.
- •Installed and operational in 2013.

Shielding

- Cylindrical shielding of 2 inches Cu + 8 inches Pb
- Nitrogen purge at 2L/min to keep radon out, as the lab radon levels are 150 Bq/m³.

Detection Region

• Energy: 10 – 900 keV

Unshielded and Shielded Spectra (Canberra Well Detector)

Energy (keV)



17/06/2015

SN

AB

15

Canberra Well Detector Background (is the detector an ultra-low counter)

Background run completed, ~150 days.

Isotope	Counts per day
²³⁸ U	1.16
²³² Th	0.51
²²⁸ Ac	0.39
²³⁵ U	0.48
⁴⁰ K	0
²¹⁰ Pb	0

- Total backgrounds at the level of ~2.5 counts / day in regions of interest. An ultra-low detector should have a background no more than 0.3 counts/day per gamma energy.
- Calibration sources approved by SNOLAB and efficiency measurements up to ~900 keV have been completed.
- Samples for DAMIC, DEAP, PICO and SNO+ have been counted or are in progress.

Canberra Well Detector Sensitivity

Isotope	Sensitivity for Standard Size Samples	Sensitivity for Standard Size Samples
²³⁸ U (↑ ²²⁶ Ra)	0.05 mBq/kg	4 ppt
²³⁸ U (↓ ²²⁶ Ra)	0.08 mBq/kg	6 ppt
²²⁸ Ac	0.2 mBq/kg	49 ppt
²³² Th	0.4 mBq/kg	98 ppt
²³⁵ U	0.02 mBq/kg	35 ppt
²¹⁰ Pb	0.15 mBq/kg	12 ppt

SN

AB

Measurements To Date For Each Experiment

Experiment or Laboratory	Total (2013 - Today)				
SNO+	30				
SNOLAB	1				
DEAP	1				
COUPP / PICO	13				
DAMIC	1				
Total	46				
Calibrations & Tests	15				
Counting time per sample averages 7 days.					

17/06/2015

SN

AB

Electrostatic Counting System



9 counters located at SNOLAB,1 on loan to LBL (EXO),1 on loan to U of A (DEAP).

Originally built for SNO, now used primarily by EXO. However, these counters are owned by SNOLAB so samples can be measured for other experiments.

Measures ²²²Rn, ²²⁴Ra and ²²⁶Ra levels.

Sensitivity Levels are:

²²²Rn: 10⁻¹⁴ gU/g

²²⁴Ra: 10⁻¹⁵ gTh/g

²²⁶Ra: 10⁻¹⁶ gU/g

Work is ongoing to improve sensitivity even further.

17/06/2015

SN

Alpha Beta Counting System



Currently located at the SNOLAB hot lab at LU so that radioactive spike sources can be measured for SNO+.

Sensitivity for 238 U and 232 Th is ~ 1 mBq assuming that the chains are in equilibrium.

17/06/2015

SN



SNOLAB maintains a database in a spreadsheet format for each experiment. https://www.snolab.ca/users/services/gamma-assay

The table shows data from the standard gamma searches: ²³⁸U, ²³⁵U, ²³²Th, ⁴⁰K ¹³⁷Cs, ⁶⁰Co.

While searching for the above gammas, we also search for any other peaks in the spectrum between 100 keV and 2800 keV, For example, ⁵⁴Mn is usually observed in steel.

SNOLAB is a member of the Assay and Acquisition of Radiopure Materials (AARM) Collaboration, which has developed the Community Material Assay Database radiopurity.org.

Future Low Background Counters and Facilities

A new dedicated space will be constructed at SNOLAB for a low background lab located in the South Drift.

This drift is isolated from other drifts and is inaccessible to large equipment. This will help reduce micro-seismic noise which can effect Ge detectors.

Increased air flow and perhaps other radon reduction techniques will be used. It is known that the compressed air from surface has substantially less radon than the lab air and can be used to reduce radon levels from 135-150 Bq/m³ to 1-5 Bq/m³.

Space can accommodate up to 5 Ge detectors, XRF, radon emanation chamber and have room for other types of counters which would benefit from low-cosmic ray background.

Engineering design drawings are now in progress.



South Drift

AB

SN



Additional Low Background Counters Coming Soon

Two additional high purity germanium detectors will be installed.

1. SNOLAB Canberra 400 cm³ coaxial detector acquired in 2011 and refurbished into an ultra-low counter in 2013 to be installed, the shielding apparatus is currently being designed.



2. Bern Coax detector, detector is currently at Bern and will be relocated to SNOLAB. This detector has been used extensively by the EXO experiment. It is expected to be shipped to SNOLAB this summer.





Future Improvements and R&D

- Improved neutron shields (detector response, spectrum)
- Improved material selection (more sensitive, better radiopurity e.g. PbWO₄ with archaeological lead)
- Active shielding

AB

SN

- Going deeper underground
- Storage of freshly made construction materials underground
- Multisegmented crystals or multiple crystals
- Collaboration with producers (e.g. depleted Ge, crystal growing, Cu electroforming underground)
- Reorganisation and optimisation of existing screening facilities is necessary, because they are costly and measurement times can be rather lengthy.
- Harmonisation of how to report data and intercomparison programs for ultra low-level measurement techniques.

Applications Beyond Particle Physics

- Ultra low-level chemistry
- Particle astrophysics (material and techniques applicable to rare events experiments)
- Space science (e.g. micro meteorites, Mars samples, cosmic activation products, comet tail samples)
- Atmospheric samples (very) short lived isotopes, radionuclide composition, stratospheric samples)
- Ocean samples (e.g. deep ocean water ⁶⁰Fe)
- In general application of low background techniques to interdisciplinary fields:
 - Low-level environmental radioactivity measurement and monitoring
 - Radiodating (extension of determined ages towards the past)
 - Geophysics (palaeoseismology, palaeogeology, sedimentation)

SN

Summary

• SNOLAB PGT HPGe low background counting system has run continuously for the past since 2005 and has counted 480 samples so far.

Counting queue is usually long.

The counter is available for all SNOLAB experiments and can be made available to non-SNOLAB experiments upon request (eg. DM-ICE, DRIFT).

• Two Canberra Ge detectors were delivered to SNOLAB, but each needed to be refurbished.

The Canberra Well detector is now in full operation and 45 samples have now been counted.

The Canberra Coax detector is underground and engineering drawings of the shielding design are in progress.

The EXO Bern detector will be shipped to SNOLAB within weeks.

- Specialized counting can be done using the Electrostatic Counters, Alpha-Beta Counters and materials can be emanated for Radon.
- New low background counting lab will be constructed at SNOLAB, final design drawings are now underway.

17/06/2015

SN

ΔB